



Productos para Laboratorio

www.acefesa.com

Agarose, Low EEO, NBS


All purpose Agarose for Molecular Biology

*** Low EEO (0.05-0.13)**

**** High Gel Strength (≥ 1200)**

***** Low Absorption of Staining Reagents**

Agarose - Low EEO is an all-purpose Agarose ideal for nucleic acid analytical/ preparative work, blotting and routine applications. This Agarose has a low EEO (electroendosmosis) allowing for shorter electrophoretic runs without compromising resolution combined with a high mechanical resistance for easier and more reliable handling.



1% Agarose - Low EEO Gel
(stained using EtBr)

Agarose - Low EEO has low absorption of staining agents and high thermal stability due to high hysteresis (difference between gelling and melting temperatures), it has excellent transparency of the gel and high visibility. The low cost of Agarose - Low EEO makes it the ideal performance Agarose for general use in electrophoresis.

Moisture $\leq 7\%$
Ash $\leq 4\%$
EEO* 0.05-0.13
Sulphate $\leq 0.1\%$
Clarity 1.5% (NTU) ≤ 3
Gel Strength 1% (g/cm²) ≥ 1200
Gel Strength 1.5% (g/cm²) ≥ 2500

Gelling Temperature 1.5% (oC) 36 ± 1.5
Melting Temperature 1.5% (oC) 88 ± 1.5
DNase/RNase None Detected
DNA Resolution $> 1000\text{bp}$ Finely Resolved
Gel Background Very Low
*EEO (electroendosmosis)

PRODUCT CODE	DESCRIPTION	PACK SIZE	P.V.P.
AGAROSE			
NBS-AG500	Agarose - Low EEO	500g	169,50€

IVA no incluido

Para pedidos inferiores a 150€ se cargarán 10€ de gastos de gestión. Portes pagados.

OFERTA VÁLIDA HASTA EL 31/12/2018

ACEFE, S.A.



Sede Central

Calle Isaac Peral, 5
Pol. Ind. Camí Ral
08850 Gavà, SPAIN
Tel: +34 936-362-002
acefesa@acefesa.com

Delegación Levante y Sur de España

Charo Ortín López
Tel.: +34 670 906 508
comercial2@acefesa.com

Delegación Nor Oeste de España

Francisco Gutiérrez Gómez
Tel.: +34 627 449 343
comercial7@acefesa.com